# pyatuomagic

Release 0.0.1

C. Zurn, R. Bechtold, A. Lawrence, L. Xiang, S. Meza, D. Soni

## **CONTENTS:**

1	Project	1
2	Block	5
3	Subject	8
4	Calculate Quality	9
5	Rate Quality	10
6	Preprocess	11
7	Perform EOG Regression	13
8	Perform Filter	14
9	Robust Principal Component Analysis	15
Рy	Python Module Index	
In	Index	

## **ONE**

## **PROJECT**

**class** Project.**Project** (name, d\_folder, file\_ext, montage, sampling\_rate, params)

Object containing all methods for creating a new project.

#### **Parameters**

- name (str) The name of the project
- **d folder** (str) The folder where the raw data is stored
- **file\_ext** (*str*) File extension
- montage (str) Montage to be used
- **sampling\_rate** (*int*) Sampling rate for the txt file
- params (dict) Preprocessing parameters for the new project

#### **Variables**

- quality\_thresholds The thresholds to rate te quality of the datasets
- ds\_rate Sampling rate to create reduced files
- rate\_cutoffs Sampling rate to the recorded data
- config Configuration file with all the project constants
- params (dir) The default parameters to be used
- **visualization\_params** (dir) The default visualisation parameters to be used
- CGV (dict) Constant Global Variables

#### create\_ratings\_structure()

Method that creates and initializes all data structures based on the data on both data folder and results folder

#### **Parameters None**

#### Returns

- block list
- processed\_list
- block\_map
- n\_processed\_subjects
- n\_processed\_files
- n block

- current
- interpolate\_list
- good\_list
- bad\_list
- ok list
- not rated list

#### static dir\_not\_hiddens(folder, extn)

Returns the list of files in the folder, excluding the hidden files

#### **Parameters**

- folder The folder in which the files are listed
- extn Extension of the raw file

Returns List of files that are not hidden

**Return type** files

#### get\_current\_block()

Returns the block pointed by the current index

**Parameters None** 

Returns Current block

Return type block

#### get\_quality\_ratings(cutoffs)

Parameters cutoffs (dict) – cutoffs on which the ratings will be decided

**Returns** ratings – list of block-wise ratings

Return type list

#### get\_rated\_count()

**Parameters None** 

Returns

Return type Count for no. of blocks that are yet to be rated

#### interpolate\_selected()

Interpolates all the channels selected to be interpolated

**Parameters None** 

Returns

Return type None

#### list\_preprocessed\_subjects()

Method that lists all folders in the results folder

**Parameters None** 

**Returns** listb – List of all folders in the results folder

Return type list

#### list\_subject\_files()

Method that lists all folders in the data folder

**Parameters None** 

Returns listb – List of all folders in the data folder

Return type list

#### static list\_subjects(root\_folder)

Returns the list of subjects in the folder

**Parameters root\_folder** (str) – The folder in which the subjects are looked for

Returns subjects – List of subjects in the root folder

Return type list

#### static make\_rating\_manually(block, q\_rate)

Returns q\_rate if the block is not rated manually.

#### **Parameters**

- **q\_rate** (*float*) The rate to be returned
- block Block for which the rate is returned

**Returns** Return q\_rate if the block is not manually rated. If it is rated manually return 'Manually rated'

Return type rate

#### preprocess\_all()

Preprocess all files in the data folder of the project

**Parameters None** 

**Returns** 

Return type None

#### save\_project()

saves the project information to a JSON file

**Parameters None** 

Returns

Return type None

#### set\_data\_folder(folder)

Sets the path folder where the data is stored

**Parameters folder** (*str*) – Path to the raw data folder

**Returns** self.data\_folder – Path to the raw data folder

Return type str

#### set\_name (name)

Sets the name of a new project

**Parameters name** (str) – Name of the current project

**Returns** self.name – Name of the new project

Return type str

#### set\_results\_folder(folder)

Sets the path folder where the results will be stored

**Parameters folder** (*str*) – Path where is the data stored

**Returns** self.results\_folder – Path to the results folder

Return type str

to\_be\_interpolated\_count()

**Parameters None** 

Returns

Return type Count for no. of blocks that are yet to be interpolated

update\_project (preprocessed)

Updates the project information dictionary with each blocks information

Parameters preprocessed – Preprocessed files

Returns

Return type None

update\_rating\_lists(block)

Updates the rating lists according to the rating of the block.

Parameters block – block for which the rating list has to be updated

Returns

Return type None

#### **TWO**

## **BLOCK**

#### class Block .Block (root\_path, data\_filename, project, subject)

Object for all operations on an individual dataset.

Initialized using the name and path of the raw data. Preprocess, interpolate, rate for quality, and store those files.

#### **Parameters**

- **root\_path** (*str*) root directory of the BIDS project
- data filename (str) BIDS filename with extension
- **project** (*object*) project object to which this block belongs
- **subject** (*object*) subject object to which this block belongs

#### Variables

- unique\_name (str) raw file name minus extension, used for saving results as well
- **file\_ext** (str) raw file extension
- params (dict) parameters for preprocessing and calculating quality metrics
- **sampling\_rate** sampling rate of raw data file
- result\_path (str) directory path to where results are stored for the block
- rate (str) current rating of the file (good, bad, ok, not rated)
- to be interpolated (list) list of channel indices that are to be interpolated
- auto bad chans (list) list of channel indices detected as bad
- final\_bad\_chans (list) list of channel indices determined to be bad after checks
- quality\_scores (dict) contains all metrics of quality calculated for the dataset
- **times\_committed** (*int*) used to track how many changes were made to the evaluation of the data

#### preprocess()

run the block through preprocessing steps, calc quality scores, save files, write log

#### interpolate()

interpolate the dataset, update quality scores and rating, save files, write log

#### find result path()

Identifies the directory path pointing to where results stored

Following BIDS requirements, we only have either the subject folder or both subject and session.

#### Parameters none

Returns result\_path - location of results files within BIDS folder

Return type str

#### interpolate()

Interpolates bad channels to create new data and updates info

Parameters none

Returns

Return type none

#### load\_data()

Load raw data from BIDS folder

Allowing for a number of extensions, loads file

Parameters none

**Returns** 

Return type raw MNE object

#### preprocess()

Preprocesses the raw data associated with this block

none

results: dict dictionary containing all the new updates to the block and the preprocessed array

#### save\_all\_files (results, fig1, fig2)

Save results dictionary and figures to results path

#### **Parameters**

- results MNE raw object with info attribute containing
- fig1 Figure of ??
- fig2 Figure of ??

#### Returns

Return type none

#### update\_rating(update)

Takes update about ratings and stores in object

From project level object, get an update on rating info.

**Parameters update** (*dict*) – dictionary of updates

Returns

Return type none

#### update\_rating\_from\_file()

Updates block information from the file currently stored

Checks for results file, if it's there, and informaation, we update. No direct returns, but updates block fields.

Parameters none

**Returns** 

## Return type none

## write\_log(updates)

Writes a log for all of the updates its making/actions performed :Parameters: updates (dict)

### Returns

Return type Updates in log file

## THREE

## **SUBJECT**

#### class Subject.Subject(data\_folder)

SUBJECT is a class representing each subject in the dataFolder. A Subject corresponds to a folder, which contains one or more Blocks. A Block represents a raw file and it's associated preprocessed file, if any (See Block).

#### static extract name(address)

extract the name of the subject from subject data folder. :Parameters: address (str) – path of the subject folder

Returns name – name of the subject

Return type str

#### result\_path(data\_folder)

finds the result folder path for the corresponding subject data folder according to the BIDS folder hierarchy. :Parameters: **data\_folder** (*str*) – subject data folder

Returns result – corresponding result folder

Return type str

#### update\_addresses (new\_data\_path, new\_project\_path)

This method is to be called to update addresses in case the project is loaded from another operating system and may have a different path to the dataFolder or resultFolder. This can happen either because the data is on a server and the path to it is different on different systems, or simply if the project is loaded from a windows to a iOS or vice versa. :Parameters: \* new\_data\_path (str) – updated path of the data folder

• **new\_project\_path** (*str*) – updated path of the project folder

Returns

Return type None

## **FOUR**

## **CALCULATE QUALITY**

Calculates four essential quality metrics of a data set and returns a dictionary of the quality metrics and thresholds/settings used

**The four quality metrics:** Overall high amplitude data points Timepoints of high variance across channels Ratio of bad channels Channels of high variance across time

Mean absolute voltage is also calculated.

#### **Parameters**

- Data (np. ndarray) a channels x timepoints data array of EEG
- bad\_chans (list) a list of the numbers of bad channels
- **overallThresh** (float) overall threshold for rejection of ...
- **timeThresh** (*float*) threshold for rejecting time segments
- **chanThresh** (float) threshold for rejecting channels
- apply\_common\_avg (bool) if average referencing should be applied

#### Returns

quality\_metrics: a dictionary of the quality metrics and thresholds/settings used The four quality metrics: Overall high amplitude data points Timepoints of high variance across channels Ratio of bad channels Channels of high variance across time

Mean absolute voltage also calculated

#### Return type dict

### RATE QUALITY

```
rateQuality.rateQuality (quality_metrics: dict, overall_Good_Cutoff: float = 0.1, overall_Bad_Cutoff: float = 0.2, time_Good_Cutoff: float = 0.1, time_Bad_Cutoff: float = 0.2, bad_Channel_Good_Cutoff: float = 0.15, bad_Channel_Bad_Cutoff: float = 0.3, channel_Good_Cutoff: float = 0.15, channel_Bad_Cutoff: float = 0.3)
```

Rates datasets, based on quality measures calculated with calcQuality().

The possible ratings: Good overall rating Regular overall rating Bad overall rating

#### **Parameters**

- quality\_metrics (dict) a dictionary containing the quality metrics to rate the dataset.
- **overall\_Good\_Cutoff** (*float*) cutoff for "Good" quality based on overall high amplitude data points [0.1].
- **overall\_Bad\_Cutoff** (*float*) cutoff for "Bad" quality based on overall high amplitude data point [0.2].
- **time\_Good\_Cutoff** (*float*) cutoff for "Good" quality based on time points of high variance across channels [0.1].
- **time\_Bad\_Cutoff** (*float*) cutoff for "Bad" quality based on time points of high variance across channels [0.2].
- bad\_Channel\_Good\_Cutoff (*float*) cutoff for "Good" quality based on ratio of bad channels [0.15].
- bad\_Channel\_Bad\_Cutoff (*float*) cutoff for "Bad" quality based on ratio of bad channels [0.3].
- **channel\_Good\_Cutoff** (*float*) cutoff for "Good" quality based on channels of high variance across time [0.15].
- **channel\_Bad\_Cutoff** (*float*) cutoff for "Bad" quality based on channels of high variance across time [0.3].

**Returns dataset\_qualification** – a dictionary indicating is the dataset if "Good" = 100, "Regular" = 50 or "Bad" = 0.

Return type dict

## **PREPROCESS**

#### class preprocess.Preprocess(eeg, params)

Preprocess class for pyautomagic preprocessing pipeline: preprocess performs pyprep's prep\_pipeline, then filters the eeg data, performs optionaleog\_regression, preforms RPCA to remove noise, and lastly can output plots of the data at various stages in the process.

#### **Parameters**

- eeg (mne.io.Raw) mne raw object containing eeg data and all the data's information.
- **params** (*dict*) dictionary of parameters <default> params = { 'line\_freqs' : 50,

```
'filter_type': 'high', 'filt_freq': None, 'filter_length': 'auto', 'eog_regression': False, 'lam': -1, 'tol': 1e-7, 'max_iter': 1000, 'interpolation_params': {'line_freqs': 50,
```

```
'ref_chs': eeg.ch_names, 'reref_chs': eeg.ch_names, 'montage': 'standard_1020'}
```

#### Variables

- eeq (mne.io.Raw) mne raw object containing eeg data and all the data's information.
- eog (mne.io.Raw) mne raw object containing eog data and all the data's information.
- bad\_chs (List) list of the names of all the detected bad channels
- params (dict) dictionary of parameters described above
- index (numpy.array) array of bad channel indices
- **filtered** (mne.io.Raw) mne raw object containing eeg data after filtering
- **eeg\_filt\_eog** (*mne.io.Raw*) mne raw object containing eeg data after filtering, and eog\_regression
- **eeg\_filt\_eog\_rpca** (*mne.io.Raw*) mne raw object containing eeg data after filtering, eog\_regression, and robust PCA (final cleaned data)
- noise (numpy.array) array of the noise removed from rpca
- automagic (dict) automagic holds information about the progress of the pipeline
- **fig1** (*matplotlib.pyplot.figure*) figure of 6 subplots at each stage of the preprocess pipeline
- fig2 (matplotlib.pyplot.figure) Figure of the final cleaned eeg data

#### perform\_prep()

Calls pyprep's PrepPipeline and detects bad channels in the data.

#### perform\_filter():

Performs initial filter (high, low, or band-pass) and removes line noise.

#### def perform\_eog\_regression

If requested, it will remove artifact from eog data.

#### perform RPCA()

Uses Robust Principal Component Analysis to remove noise from the data.

#### plot(self, show=True):

Outputs plots of data at each point in the preprocess pipeline.

#### fit(self):

Perform the full preprocessing pipeline for pyautomagic (modeled from matlab's automagic package).

#### fit()

Fit Perform the full preprocessing pipeline for pyautomagic (modeled from matlab's automagic package).

#### Returns

- **eeg\_filt\_eog\_rpca** (*mne.io.Raw*) Corrected Data
- fig1 (matplotlib.pyplot.figure) figure of 6 subplots at each stage of the preprocess pipeline
- fig2 (matplotlib.pyplot.figure) Figure of the final cleaned eeg data

#### perform\_RPCA()

Uses Robust Principal Component Analysis to remove noise from the data.

#### Returns

- **eeg\_filt\_eog\_rpca** (*numpy.array*) eeg data after filtering, eog\_regression, and robust PCA (final cleaned data)
- noise (numpy.array) array of the noise removed from rpca

#### perform\_eog\_regression()

If requested, it will remove artifact from eog data.

Returns eeg\_filt\_eog - Filtered eeg data, with eog regression

Return type mne.io.Raw

#### perform\_filter()

Performs initial filter (high, low, or band-pass) and removes line noise.

Returns filtered - Filtered eeg data

Return type mne.io.Raw

#### perform prep()

Calls pyprep's PrepPipeline and detects bad channels in the data.

Returns bad\_chs - List of all the bad channel names

Return type List

#### plot (show=True)

Outputs plots of data at each point in the preprocess pipeline.

#### Returns

- fig1 (matplotlib.pyplot.figure) figure of 6 subplots at each stage of the preprocess pipeline
- fig2 (matplotlib.pyplot.figure) Figure of the final cleaned eeg data

**SEVEN** 

## PERFORM EOG REGRESSION

 ${\tt perform\_EOG\_regression.perform\_EOG\_regression} \ (\textit{EEG}, EOG)$ 

The artifacts due to EOG activity are removed from the EEG data using the subtraction method that relies on the linear transformation of the EEG signal

#### **Parameters**

- **EEG** (*np.ndarray*) The input EEG data
- EOG (np.ndarray) The input EOG data

Returns clean\_EEG - Cleaned EEG signal from EOG artifacts

Return type np.ndarray

#### References

[1] Pedroni, A., Bahreini, A., & Langer, N. (2019). Automagic: Standardized preprocessing of big EEG data. Neuroimage, 200, 460-473. doi: 10.1016/j.neuroimage.2019.06.046

## **EIGHT**

## **PERFORM FILTER**

performFilter.performFilter(*EEG*, *sfreq*, *filter\_type=None*, *filt\_freq=None*, *filter\_length='auto'*)

This function filters EEG data using Hamming windowed sinc FIR filter with input filter type and parameters.

#### **Parameters**

- **EEG** (ndarray, shape (...,  $n\_times$ )) Input EEG data to be filtered.
- **sfreq** (*float* | *None*) The sample frequency in Hz.
- filter\_type (str) The filter type, can only take 'low', 'high' or 'notch', defaults to None
- filt\_freq (float | None) The filter frequency in Hz, defaults to None
- **filter\_length** (*str* | *int*) Length of the FIR filter to use, defaults to 'auto'.

Returns The filtered EEG data.

**Return type** ndarray, shape (..., n\_times)

## ROBUST PRINCIPAL COMPONENT ANALYSIS

rpca.**rpca**(*M*, *lam=-1*, *tol=1e-07*, *maxIter=1000*)

Perform Robust Principle Component Analysis:

Performs a Robust Principal Component Analysis on the EEG data with the specified parameters: Lamda, Tolerance, and Maximum number of Iterations. The function outputs the EEG data with the noise removed as well as the nosie that was removed.

Adapted from Cian Scannell - Oct-2017 (https://github.com/cianmscannell/RPCA) Computes rpca separation of M into L(low rank) and S(Sparse) using the parameter lam this uses the alternating directions augmented method of multipliers.

#### **Parameters**

- M (*npumpy.darray*) 1st parameter, EEG Data (must include)
- lam (double) 2nd parameter, Lamda paramter for RPCA (default = 1/(sqrt(# of Colunms))
- tol (double) 3rd parameter, Tolerance (defalut = 1e-7) RPCA param
- maxIter (*int*) fourth parameter, Maximum Iterations (deafult = 1000)

#### Returns

- **Data** (*npumpy.darray*) Corrected Data (Low rank matrix)
- Error (npumpy.darray) Noise removed from the data (Sparse Matrix) note: M = L + S

#### rpca.soft thres (x, eps)

Cian Scannell - Oct-2017 Soft thresholds a matrix x at the eps level i.e  $ST(x,eps)_{ij} = sgn(x_{ij}) max(|x_{ij}| - eps, 0)$ 

#### **Parameters**

- x (npumpy.darray) first parameter, values to be thersholded
- eps (double) second parameter, thershold

**Returns np.multiply(a,b)** – thersholded values, where anythign under the threshold was set to zero **Return type** npumpy.darray

## **PYTHON MODULE INDEX**

```
b
Block,5

C
calcQuality,9

p
perform_EOG_regression,13
performFilter,14
preprocess,11
Project,1

r
rateQuality,10
rpca,15

S
Subject,8
```

## **INDEX**

В	M
Block (class in Block), 5 Block (module), 5	<pre>make_rating_manually() (Project.Project static</pre>
С	P
<pre>calcQuality (module), 9 calcQuality() (in module calcQuality), 9 create_ratings_structure() (Project.Project</pre>	<pre>perform_EOG_regression (module), 13 perform_EOG_regression() (in module per- form_EOG_regression), 13 perform_eog_regression() (prepro- cess.Preprocess method), 12 perform_filter() (preprocess.Preprocess method),</pre>
<pre>dir_not_hiddens() (Project.Project static method), 2</pre>	perform_prep() (preprocess.Preprocess method), 11, 12
<pre>E extract_name() (Subject.Subject static method), 8 F find_result_path() (Block.Block method), 5</pre>	perform_RPCA() (preprocess.Preprocess method), 12 performFilter (module), 14 performFilter() (in module performFilter), 14 plot() (preprocess.Preprocess method), 12 Preprocess (class in preprocess), 11
Git () (preprocess.Preprocess method), 12  Get_current_block() (Project.Project method), 2  get_quality_ratings() (Project.Project method),	preprocess (module), 11 preprocess() (Block.Block method), 5, 6 preprocess_all() (Project.Project method), 3 Project (class in Project), 1 Project (module), 1
get_rated_count() (Project.Project method), 2  interpolate() (Block.Block method), 5, 6 interpolate_selected() (Project.Project method), 2	R rateQuality (module), 10 rateQuality() (in module rateQuality), 10 result_path() (Subject.Subject method), 8 rpca (module), 15 rpca() (in module rpca), 15
L  list_preprocessed_subjects()	S save_all_files() (Block.Block method), 6 save_project() (Project.Project method), 3 set_data_folder() (Project.Project method), 3 set_name() (Project.Project method), 3 set_results_folder() (Project.Project method),

## Τ

```
to_be_interpolated_count() (Project.Project method), 4
```

## U

```
update_addresses() (Subject.Subject method), 8
update_project() (Project.Project method), 4
update_rating() (Block.Block method), 6
update_rating_from_file() (Block.Block method), 6
update_rating_lists() (Project.Project method),
```

## W

write\_log() (Block.Block method), 7

Index 18